

# NEUTONE

## Digital Display Muscle(Soft Tissue) Hardness Tester



### Special Features

Result of using sophisticated one-chip microcomputer, The NEUTONE TDM-Z1(RB) is the faster, easier, smarter way to measure Muscle(Soft Tissue) hardness.

The NEUTONE TDM-Z1(RB) is operated by rechargeable battery (no need exchange battery like Dry-cell operation).

The NEUTONE TDM-Z1(RB) features automatically display average measuring value (2-9 times), stores measuring value, also exports to your personal computer as CSV file of measuring value.

For measuring the Muscle(Soft Tissue) hardness, the spring-loaded presser is pressed downwards until trigger switch on. This ensures a constant pressure and prevents measuring errors.

- SIMPLE :** Operation is as simple as pressing down, nothing more.
- ACUURATE:** A constant-pressure uniform pressure preventing false reading due to differences between operators.
- REASONABLE:** Very affordable price but no compromised quality.



28 Sept. 2010

**TRY-ALL**

Try-All Corporation Phone 81-43-270-0588 Fax 81-43-277-3447  
Addr. 2-15-5-901 Masago Mihama-Ku, Chiba 261-0011 Japan